

Tokyo Metropolitan Bokuto Hospital

Fujitsu Hospital Information System (HIS)

Dept/ Ward Client



LS 1



LS 2



LIMS



LWS



Control Client

CLINIEELServer

Zonal Verification
Method LES



LCS

Router
TCP/IP



Ethernet 100Mbps

VPN Line
Remote Maintenance
A&T

Central Blood
Collection



BC-ROBO787
HARN570

Sample Reception

Send out




BC-ROBO-580

Network Printers



FAX Report



For Int./Ext.
Comm

Clinical Chem/
Immunology



A&T TLA System
CLINILOG Ver.3
Unloading Mod.x 1
Loading Mod.x 1
Aliquot Mod.x 2
Decapper x 1
Coagrex800 x 2
ARCHITECT x 2
BM6050 x 3
EA07U x 2

Lumipulse f

AES320
RP System (G8)
RP System (GA08)

Network Printer



Hematology
CLINILAN HM



SE- α
SE-9000

KX21

ESR6000

Hemogram Counter



Network Color Laser Printer



Network Printer

Emergency (Urine)



Clinitek 500
UF-50

Urinalysis



UA-ROBO600
OM6030

Clinitek ATLAS
HITACHI 6800

Urine Sediment Counter



Network Printer



Statistics Management
(Supervisor)



Conference



Client for LM



Network Printer

Microbiology
CLINILAN MB



Walkaway96
Bact/Alert

COBAS Amplicor
MGIT

Network Color Laser Printer



Physiology



Radiometer Blood Gas
ABL505


Transfusion Management
CLINILAN BT




AutoVueInnova



Scanner



Network Color Laser
Printer



Network Printer

NICU



Bayer 860