

Goto Central Hospital

NEC Ordering System

DB Server (Main)

DB Server (Sub)

Transaction Server (Main)


Transaction Server (Sub)

CLINIEEL Server

Communication Server

On-line Controller


Central Blood Collection Room



Test Tube Preparation System
BC-ROBO585 (For 8 Types)

Horn Cup Labeler HARN-570


Clinical Chemistry



Laboratory Automation System
CLINILOG Ver.2.0

Laboratory Automation System
1. Immunology Analyzer (ARCHITECT)
2. Clinical Chemistry Analyzer (BM-8)
3. Clinical Chemistry Analyzer (BM-1250)

Reception




IP100


SR901

OMR


Hematology



Blood Cell Counter
KX-21




Hematology Sys. Sysmex SIS




Urinalysis Sub Sys. Sysmex SIS

Transfusion




BL/Power On-line

Occult Blood Analyzer
JIA-HB2010



HbA1c Analyzer G7



Blood Sedimentation
Quick 15

Hematology Sub System

Blood Cell Counter
Smear Creator XE2100 a

Coagulation Analyzer
CA-1500

WBC Differentiation Sysmex

WBC Differentiation Sysmex


Urinalysis Sub System

Urine Analyzer
Clinitek Atlas & 500

Sedimentation Analyzer
UF-100 Sysmex

Sedimentation Differentiation
Sysmex


Serum




Immunology AXSYM

Protein Fraction CTE780


for Reference



LP9300



FD Preparation for Outsourcing



Modem

A&T

ISDN Remote Maintenance

Ethernet 100Mbps