

# Tokyo Metropolitan Komagome Hospital

**Hospital Information System (HIS)**

Dept/ Ward Client

**Main Server** **Sub Server** **Main Server 2** **Image Control** **WEB** **Expert Engine**

**Control Client** **Communication Server** **Router TCP/IP**

Ethernet 100Mbps

Dedicated line Remote Maintenance A&I

**Central Blood Collection/ Reception**

OMR

HRN570

BC-ROBO787

**System Monitoring**

**Network Printers**

**Immunology**

BN II

Lumipulsef

IMX

**Clinical Chemistry/ Emergency**

OMR

CLINILOG Ver.2 TLA System

250BloodGasSystem

250BloodGasSystem

MAOS

CTE-1800

CYTOMICS FC500

**FAX**

In-hospital line

**Send out**

**Hematology CLINILAN /HM**

XE-2000

SE-9000 α

SP-100

R-3500

ESR6000

**Hemogram Counter**

Color Laser Printer

**Urinalysis**

ATLAS5001X

**Urine Sediment Counter**

**Network Printers**

**Statistics Management**

**Technologists**

Laser Printer

**Operation**

Blood Gas Analyzer 860

**Microbiology/ Infection Control CLINILAN /MB**

Walk Away

BacTAlert

TDX

TDX

Cobas Amplicor

Cobas Amplicor


MB/Bact240

MB/Bact240

Wellreader SK601

**Color Laser Printer**

# 2F Laboratory (Server room installed)



Main Server

Sub Server

Ethernet 100Mbps

## Existing Laboratory System CLINILAN Series



TCP/IP

## Fujitsu Hospital Information System (HIS)

Dept/ Ward Client




### Blood Transfusion

#### Reception



Barcode Printer

#### Transfusion



Scanner (A4)



#### Autologous Blood Transfusion




Color Inkjet Printer


#### Statistics Management



#### Blood Product Management




#### Color Laser Network Printers



#### Blood Product Ordering



MODEM



FAX for Blood Product by Japanese Red Cross Society